


OPTIMUS

Most optimized for clinician



OPTIMUS MINI SCREW



Most optimized Mini Screw system for clinicians



OSTEONIC Co., Ltd.

505-3, Ace Techno Tower 3rd, Guro-dong, Guro-gu, Seoul, Korea
TEL: +822 6052 8886 FAX: +822 6326 7001 www.osteonic.com

ORT - E - 1208




Optimus Mini Screw system is most optimized for clinicians.


It provides clinicians with fast and effective tools to treat a full range of orthodontic cases. Optimus Mini Screws can be placed in the Maxilla and the Mandible area with purpose of offering a true alternative to conventional anchorage methods.


- Self-tapping and self-drilling to minimize the need for tissue punches or pilot drills
- Tissue-suppression stops to prevent tissue overgrowth
- Smooth and polished transmucosal collars prevent infection
- Buttress thread design enhances pull out strength
- Rounded Hexagonal head design helps to smoothly remove the shaft from screw head
- Dual cutting threads remove bone debris in high density bone
- Flute thread helps an initial drill with removing bone debris in high density bone

OPTIMUS BLOCK & INSTRUMENTS

Part	Order No.	Description
	TIN-006	Sterilization & Logistic Block
	IOC-010	Driver Handle
	DDS-L01	Driver Shaft, Double & Wing Head Type
	DDS-L02	Driver Shaft for Handle, Button Head Type
	DDS-S01 DDS-S02	Driver Shaft for Handpieces, Double & Wing Head Type
	DDS-S03 DDS-S04	Driver Shaft for Handpieces, Button Head Type

Optimus Screws

Double Head Type			
Double Head Type	Order No.	Dia.	L(mm)
	D11406	1.4	6
	D11407	1.4	7
	D11408	1.4	8
	D11606	1.6	6
	D11607	1.6	7
	D11608	1.6	8
	D11609	1.6	9
	D11806	1.8	6
	D11807	1.8	7
	D11808	1.8	8
	D11809	1.8	9
	D11810	1.8	10
	D12007	2.0	7
	D12008	2.0	8
	D12009	2.0	9
	D12010	2.0	10
	D12011	2.0	11
	D12012	2.0	12

Button Head Type			
Button Head Type	Order No.	Dia.	L(mm)
	B11406	1.4	6
	B11407	1.4	7
	B11408	1.4	8
	B11606	1.6	6
	B11607	1.6	7
	B11608	1.6	8
	B11609	1.6	9
	B11806	1.8	6
	B11807	1.8	7
	B11808	1.8	8
	B11809	1.8	9
	B11810	1.8	10
	B12007	2.0	7
	B12008	2.0	8
	B12009	2.0	9
	B12010	2.0	10
	B12011	2.0	11
	B12012	2.0	12

Wing Head Type			
Wing Head Type	Order No.	Dia.	L(mm)
	W11406	1.4	6
	W11407	1.4	7
	W11408	1.4	8
	W11606	1.6	6
	W11607	1.6	7
	W11608	1.6	8
	W11609	1.6	9
	W11806	1.8	6
	W11807	1.8	7
	W11808	1.8	8
	W11809	1.8	9
	W11810	1.8	10
	W12007	2.0	7
	W12008	2.0	8
	W12009	2.0	9
	W12010	2.0	10
	W12011	2.0	11
	W12012	2.0	12



- Non-Sterilized Optimus Mini Screw : Just use the Order-No.
- Sterilized Optimus Mini Screw : Add "S" to Order No Ex: W11406S